

5x DNA Loading Buffer Red



Cat. No.: A608104

A608104

Cat. No.	Product	ID No.	Cap colour	Volume
A608104	DNA Loading Buffer Red	5400100	Red	5 x 1 ml

Features

- Suitable for agarose and SDS DNA gels
- Suitable for TAE, TBE, SB and LB electrophoresis buffers
- 5x ready-to-use formulation
- Ficoll-based for convenient storage on RT

General Description

DNA loading buffers are used to load DNA samples to agarose or SDS DNA gels for gel electrophoresis.

DNA loading buffers serve three main purposes:

Firstly, they add density to the DNA samples, so the DNA sinks down into the well instead of floating up and mixing with the running buffer. To achieve this, high density reagents like glycerol, sucrose or Ficoll are added.

Secondly, loading buffers add visibility to the DNA sample. Loading buffers contain a coloured tracking dye to allow control of proper DNA sample loading.

And thirdly, loading buffers provide one or more tracking dyes to monitor the progress of DNA migration on the gel (figure 1). Commonly used tracking dyes are xylene cyanol FF, cresol red, bromophenol blue, and orange G, each migrating at a characteristic size (table 1).

Ampliqon offers 4 different loading buffers to allow the customer to choose the optimal system for a specific task. The loading buffers are formulated as a 5x solutions containing Ficoll, Tris-buffer, EDTA and either Xylene cyanol FF, Cresol Red, Bromophenol Blue or Orange G as tracking dye. For a 10 µl loading volume, add 2 µl 5x Loading Dye to 8 µl of your DNA sample, mix well and load on a gel.

Kit Components

5x Loading Buffer Red

15 % Ficoll 400, 10 mM Tris-HCl pH 8.0, 50 mM EDTA, 0.05 % Cresol Red.

Recommended Storage and Stability

Long term storage at -20 °C. Product expiry at -20 °C is stated on the label.

Option: Store at +4 °C for up to 6 months.

Quality Control

Each lot of Loading Buffer is functionally tested on an agarose gel.

How to choose the right loading buffer

Tracking dyes

The choice of the tracking dye is dependent on the size of the DNA fragments one wants to run. In general, the front of the tracking dye should not run at the size of the DNA fragments because especially dark dyes obscure the DNA bands.

It is best to choose a tracking dye that runs in front of the DNA fragments to analyze.




Figure 1: The three Ampliqon Loading Buffers were run on a 1 % agarose gel along with a DNA marker. The gel was pictured on white light and UV light and the dyes were correlated to the marker.

Table 1: Position of dye front of some tracking dyes on a 1 % agarose gel:

Dye:	Front at*:	Color:
Xylene cyanol FF	4000 – 5000 bp	Cyan
Cresol red	1000 – 2000 bp	Red
Bromphenol blue	300 – 500 bp	Blue
Orange G	50 – 150 bp	Orange

* The position of the dyes is dependent on the type of agarose, the percentage of the gel and the buffer used.

Related Accessories

Accessories	Cat. No.
PCR Grade Water, 6 x 5 ml	A360056
5x DNA Loading Buffer Red, 5 x 1 ml	A608104
5x DNA Loading Buffer Blue, 5 x 1 ml	A608204
5x DNA Loading Buffer Orange, 5 x 1 ml	A608304
High Range DNA Ladder, 200-12000 bp 100 lanes	A610141
Low Range DNA Ladder, 100-1000 bp 100 lanes	A610241
PCR DNA Ladder, 100-3000 bp 100 lanes	A610341

Related Products

Taq Polymerase (500 units) *	Cat. No.
Taq DNA Polymerase 5 U/μl	A110003
• with 10x Ammonium Buffer	A111103
• 5x PCR Buffer RED	A111803
Taq DNA Polymerase 5 U/μl, RED	A200003
• with 10x Ammonium Buffer	A201103
Taq DNA Polymerase 5 U/μl, glycerol free	A100003
• with 10x Ammonium Buffer	A101103

Hot Start Polymerase (500 units) *	Cat. No.
TEMPase Hot Start DNA Polymerase, 5 U/μl	A220003
• with 10x Ammonium Buffer	A221103
• 5x PCR Buffer RED	A221803
TEMPase Hot Start DNA Polymerase, glycerol free 5 U/μl	A240003
• with 10x Ammonium Buffer	A241103

High Fidelity - Proof reading (500 units) **	Cat. No.
AccuPOL DNA Polymerase 2.5 U/μl	A210003
• with 10x Ammonium Buffer	A211103

*Available in kits including one or two buffers (Ammonium Buffer, Standard Buffer or Combination Buffer). **AccuPOL only available in kits with Ammonium Buffer. All kits include extra 25 mM MgCl₂.

Buffers for DNA polymerases *	Cat. No.
10x Ammonium Buffer, 3 x 1.5 ml	A301103
10x Standard Buffer, 3 x 1.5 ml	A302103
10x Combination Buffer, 3 x 1.5 ml	A303103
5x PCR Buffer RED, 6 x 1,5 ml **	A301810

*Ammonium Buffer, Standard Buffer and Combination Buffer are also available as Mg²⁺ free buffers, detergent free buffers and Mg²⁺ and detergent free buffers.

**For direct gel loading and visualisation.

Taq Master Mixes (500 x 50 μl reactions) *	Cat. No.
2x Master Mix, 1.5 mM MgCl ₂ final concentration	A140303
2x Master Mix RED, 1.5 mM MgCl ₂ final concentration	A180303

TEMPase Hot Start Master Mixes (500 x 50 μl reactions) *	Cat. No.
2x Master Mix A**, 1.5 mM MgCl ₂ final concentration	A230303
2x Master Mix A**BLUE, 1.5 mM MgCl ₂ final concentration	A290403

*Master mixes available also in 1.1x variants as well as 2 mM MgCl₂ variants, **Mix A is Ammonium Buffer based, also available as Mix C based on Combination Buffer.

Special Master Mixes (500 x 50 μl reactions)	Cat. No.
Multiplex 2x Master Mix, 3 mM MgCl ₂ final concentration	A260303
GC TEMPase 2x Master Mix I – for GC-rich templates	A331703
GC TEMPase 2x Master Mix II – for GC-rich templates	A332703

Real-time PCR Master Mixes (400 x 25 μl reactions)	Cat. No.
RealQ Plus 2x Master Mix for probe,	
• without ROX™	A313402
• with low ROX™	A314402
• with high ROX™	A315402
RealQ Plus 2x Master Mix Green	
• without ROX™	A323402
• with low ROX™	A324402
• with high ROX™	A325402

Ultrapure dNTPs*	Cat. No.
dNTP Mix 40 mM (2 x 500 μl): 10 mM each dA, dC, dG, dT	A502004
dNTP Set, 100 mM each: 250 μl of each dA, dC, dG and dT	A511104

*Other concentrations and Single dNTPs are available.

Loading Buffers and Ladders	Cat. No.
5x Loading Buffer Red *, 5 x 1 ml	A608104
PCR DNA Ladder **, 100 – 3000 bp, 1 x 0.5 ml	A610341

* Also available with Blue, Orange or Cyan. ** Available in different size ranges.

Reagents for *in vitro* laboratory use only.

Other product sizes, combinations and customized solutions are available. Please look at www.ampliqon.com or ask for our complete product list for PCR Enzymes. For customized solutions please contact us.

Made in Denmark
Issued 03/2017