



## rSAP

Concentration: 1 unit/ $\mu$ l

Cat. No.: A650101  
150 Units



-	rSAP 1 U/ $\mu$ l
ID No.	7200300
Cap colour	Yellow
Content	1 x 0.15 ml

MADE IN DENMARK

### Key Features

- Dephosphorylation of 5' and 3' ends of DNA and RNA
- Degradation of dNTPs prior to sequencing or SNP analysis
- Dephosphorylation prevents ligation of linearized DNA
- Dephosphorylation of DNA prior to end-labelling using T4 Polynucleotide Kinase
- Specific activity  $\geq$  2 000 U/mg

Recombinant Shrimp Alkaline Phosphatase (rSAP) is a heat labile alkaline phosphatase purified from *Pichia pastoris*, carrying the shrimp alkaline phosphatase gene from *Pandalus borealis*. rSAP is completely and irreversibly inactivated by heating at 65°C for 5 minutes, thereby making removal of rSAP prior to ligation or end-labeling unnecessary. Optimum working range for rSAP is between pH 7-9. rSAP is active in most restriction and PCR buffers. Mg<sup>2+</sup> (>1 mM) is required for activity.

### Components

#### Ampliqon rSAP in storage buffer

- rSAP 1 unit/ $\mu$ l, 25 mM Tris-HCl, pH 7.6 at 4 °C, 5 mM MgCl<sub>2</sub>, glycerol 50% (v/v)

#### Recommended Storage and Stability

Long term storage at -20 °C. Product expiry at -20 °C is stated on the label.

#### Unit Definition

One unit is defined as the amount of enzymes converting 1  $\mu$ mol of p-nitrophenyl phosphate per minute to nitrophenol and phosphate at 37°C and pH 10.4 in 0.1 M glycine buffer, 1 mM each of ZnCl<sub>2</sub> and MgCl<sub>2</sub> and 6 mM 4-nitrophenyl phosphate.

#### Quality Control

rSAP is tested for double stranded and single stranded endonuclease activity. Furthermore, the activity of rSAP is measured. rSAP is functionally tested (together with HL-ExoI) by spiking a PCR product with dNTPs and primers followed by Sanger sequencing.

#### Applications

- Removal of phosphorylated ends of DNA and RNA for downstream applications such as cloning or end-labeling of probes. Dephosphorylation prevents ligation of linearized DNA. rSAP acts on 3' protruding, 5' protruding, 5' recessed and blunt ends.
- Degradation of unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing or SNP analysis.

### Related Products

<b>rSAP 1 units/<math>\mu</math>l</b>	<b>Cat. No.</b>
• 750 Units	<b>A650103</b>
<b>HL-ExoI 20 units/<math>\mu</math>l</b>	<b>Cat. No.</b>
• 1000 Units	<b>A640103</b>
• 5000 Units	<b>A640107</b>
<b>PureIT ExoZAP PCR CleanUp Kit</b>	<b>Cat. No.</b>
• 500 reactions	<b>A630203</b>
• 5000 reactions	<b>A630207</b>
<b>PureIT ExoZAP PCR CleanUp</b>	<b>Cat. No.</b>
• 500 reactions	<b>A620603</b>
• 5000 reactions	<b>A620607</b>
<b>Taq Polymerase (500 units) *</b>	<b>Cat. No.</b>
Taq DNA Polymerase 5 U/ $\mu$ l	<b>A110003</b>
• with 10x Ammonium Buffer	<b>A111103</b>
• with 5x PCR Buffer RED	<b>A111803</b>
Taq DNA Polymerase 5 U/ $\mu$ l, RED	<b>A200003</b>
• with 10x Ammonium Buffer	<b>A201103</b>
Taq DNA Polymerase 5 U/ $\mu$ l, glycerol free	<b>A100003</b>
• with 10x Ammonium Buffer	<b>A101103</b>
<b>Hot Start Polymerase (500 units) *</b>	<b>Cat. No.</b>
TEMPase Hot Start DNA Polymerase, 5 U/ $\mu$ l	<b>A220003</b>
• with 10x Ammonium Buffer	<b>A221103</b>
• 5x PCR Buffer RED	<b>A221803</b>
TEMPase Hot Start DNA Polymerase, glycerol free 5 U/ $\mu$ l	<b>A240003</b>
• with 10x Ammonium Buffer	<b>A241103</b>
<small>*Available in kits including one or two buffers (Ammonium Buffer, Standard Buffer or Combination Buffer).</small>	
<b>Buffers for DNA polymerases *</b>	<b>Cat. No.</b>
10x Ammonium Buffer, 3 x 1.5 ml	<b>A301103</b>
10x Standard Buffer, 3 x 1.5 ml	<b>A302103</b>
10x Combination Buffer, 3 x 1.5 ml	<b>A303103</b>
5x PCR Buffer RED, 3 x 1,5 ml **	<b>A301803</b>
<small>*Ammonium Buffer, Standard Buffer and Combination Buffer are also available as Mg<sup>2+</sup> free buffers, detergent free buffers and Mg<sup>2+</sup> and detergent free buffers. **For direct gel loading and visualisation.</small>	
<b>Taq Master Mixes (500 x 50 <math>\mu</math>l reactions) *</b>	<b>Cat. No.</b>
2x Master Mix, 1.5 mM MgCl <sub>2</sub> final concentration	<b>A140303</b>
2x OptiMix CLEAR, 1.5 mM MgCl <sub>2</sub> final concentration	<b>A370503</b>
2x Master Mix RED, 1.5 mM MgCl <sub>2</sub> final concentration	<b>A180303</b>
<b>TEMPase Hot Start Master Mixes (500 x 50 <math>\mu</math>l reactions) *</b>	<b>Cat. No.</b>
2x Master Mix A**, 1.5 mM MgCl <sub>2</sub> final concentration	<b>A230303</b>
2x Master Mix A**BLUE, 1.5 mM MgCl <sub>2</sub> final concentration	<b>A290403</b>
<small>*Master mixes available also in 1.1x variants as well as 2 mM MgCl<sub>2</sub> variants, **Mix A is Ammonium Buffer based, also available as Mix C based on Combination Buffer.</small>	
<b>Special Master Mixes (500 x 50 <math>\mu</math>l reactions)</b>	<b>Cat. No.</b>
Multiplex 2x Master Mix, 3 mM MgCl <sub>2</sub> final concentration	<b>A260303</b>
GC TEMPase 2x Master Mix I – for GC-rich templates	<b>A331703</b>
GC TEMPase 2x Master Mix II – for GC-rich templates	<b>A332703</b>
<b>Real-time PCR Master Mixes (400 x 25 <math>\mu</math>l reactions)</b>	<b>Cat. No.</b>
RealQ Plus 2x Master Mix for probe,	
• without ROX <sup>TM</sup>	<b>A313402</b>
• with low ROX <sup>TM</sup>	<b>A314402</b>
• with high ROX <sup>TM</sup>	<b>A315402</b>
RealQ Plus 2x Master Mix Green	
• without ROX <sup>TM</sup>	<b>A323402</b>
• with low ROX <sup>TM</sup>	<b>A324402</b>
• with high ROX <sup>TM</sup>	<b>A325402</b>

Reagents for *in vitro* laboratory use only.

Other product sizes, combinations and customized solutions are available. Please look at [www.ampliqon.com](http://www.ampliqon.com) or ask for our complete product list for PCR Enzymes. For customized solutions please contact us.