

## TROUBLESHOOTING

The table below covers most problems and solutions hereof, experienced when handling PCR experiments.

To avoid spending a lot of time on optimisation of PCR setup we recommend the usage of Ampliqon Ammonium Buffer for most PCR applications. Ammonium Buffer is a very robust 10x PCR buffer, resulting in high yield of PCR products and minimises the need for optimisation of Mg<sup>2+</sup> and/or annealing temperatures.

Observed problem	Possible cause	Solution
<b>PCR product does not have the correct size</b>	Contamination by nucleases	<ul style="list-style-type: none"> <li>Try again with fresh reagents</li> </ul>
	Mispriming	<ul style="list-style-type: none"> <li>Test that primers do not have additional complementary regions within the template DNA</li> </ul>
	Non optimal MgCl <sub>2</sub> concentration	<ul style="list-style-type: none"> <li>Adjust MgCl<sub>2</sub> concentration as advised in product data sheet.</li> </ul>
	Non optimal annealing temperature	<ul style="list-style-type: none"> <li>Retest T<sub>m</sub> values of primers</li> </ul>
<b>Absence of PCR product</b>	Low primer specificity	<ul style="list-style-type: none"> <li>Verify that primers are complementary to the correct target sequence</li> </ul>
	Too low primer concentration	<ul style="list-style-type: none"> <li>Adjust in the range 0.1 – 1 μM</li> </ul>
	Suboptimal reaction conditions	<ul style="list-style-type: none"> <li>Optimise annealing by running a temperature gradient</li> <li>Adjust MgCl<sub>2</sub> concentration as advised in product data sheet</li> </ul>
	Poor template quality	<ul style="list-style-type: none"> <li>Test DNA using gel electrophoresis before and after addition of MgCl<sub>2</sub>.</li> <li>Check 260/280 ratio of DNA template</li> </ul>
	Missing a reaction component	<ul style="list-style-type: none"> <li>Make a new PCR mix</li> </ul>
	Inhibitors in the reaction	<ul style="list-style-type: none"> <li>Ensure that template DNA is purified or decrease sample volume.</li> </ul>
	PCR run is non optimal	<ul style="list-style-type: none"> <li>Add more cycles</li> <li>Recheck the PCR program</li> <li>Recalibrate heating block</li> </ul>
	Your template or target is complex	<ul style="list-style-type: none"> <li>For GC-rich sequences or other complex DNA targets optimize conditions using GC-rich Target kit.</li> </ul>
<b>Smears or multiple band on the gel</b>	Premature replication	<ul style="list-style-type: none"> <li>Use TEMPase Hot Start DNA Polymerase instead</li> <li>Set PCR reaction up on ice.</li> </ul>
	Too low annealing temperature	<ul style="list-style-type: none"> <li>Increase annealing temperature</li> <li>If not already using Ammonium Buffer, then shift to this buffer.</li> </ul>
	Excess primers	<ul style="list-style-type: none"> <li>Adjust in the range 0.1 – 1 μM</li> </ul>
	Non optimal MgCl <sub>2</sub> concentration	<ul style="list-style-type: none"> <li>Adjust MgCl<sub>2</sub> concentration as advised in product data sheet</li> </ul>
	Non optimal primer design	<ul style="list-style-type: none"> <li>Ensure that primers are non-complementary</li> <li>Increase length of primers</li> <li>Avoid GC-rich 3' ends</li> </ul>
	Contamination with non-template DNA	<ul style="list-style-type: none"> <li>Always use filter tips, PCR grade water.</li> <li>Use separate areas for PCR reaction setup, DNA preparation, PCR thermal cycling and gel electrophoresis</li> </ul>
	Incorrect template concentration	<ul style="list-style-type: none"> <li>Adjust template concentration as advised in product data sheet.</li> </ul>
<b>Sequence errors</b>	Low fidelity polymerase	<ul style="list-style-type: none"> <li>Use AccuPOL DNA Polymerase with Ammonium Buffer</li> </ul>
	Template DNA has been damaged	<ul style="list-style-type: none"> <li>Prepare a new DNA template</li> <li>Limit the exposure of template DNA to UV</li> <li>Lower initial heating time</li> </ul>
	Suboptimal reaction conditions	<ul style="list-style-type: none"> <li>Decrease extension time</li> <li>Decrease MgCl<sub>2</sub> concentration</li> <li>Lower the amount of cycles</li> </ul>
	Problems with nucleotide composition	<ul style="list-style-type: none"> <li>Make a fresh solution of nucleotide mix</li> </ul>