

REAL-TIME PCR MASTER MIXES




Real-time PCR Master Mixes with high efficiency and great prices


- All-in-one optimized real-time PCR Master Mixes.
- High efficiency for accurate quantification experiments.
- Convenient Hot Start capability for room temperature setup
- High specificity and stability.

RealQ Plus 2X Master Mixes from Ampliqon are easy to use all-in-one Master Mixes for both dye- and probe based detection. The mixes are available in three different levels of ROX™ reference dye and have been optimized by the Ampliqon R&D team to have high efficiency and specificity for more reliable results.


Low Standard Deviation



High Efficiency



High Specificity



	Size	Content	Cat #
	Reaction size: 25 µl		
RealQ Plus 2x Master Mix Probe High ROX™	400 Reactions	4 x 1.25 ml	A315402
	4000 Reactions	40 x 1.25 ml	A315406
RealQ Plus 2x Master Mix Probe Low ROX™	400 Reactions	4 x 1.25 ml	A314402
	4000 Reactions	40 x 1.25 ml	A314406
RealQ Plus 2x Master Mix Probe Without ROX™	400 Reactions	4 x 1.25 ml	A313402
	4000 Reactions	40 x 1.25 ml	A313406
RealQ Plus 2x Master Mix Green High ROX™	400 Reactions	4 x 1.25 ml	A325402
	4000 Reactions	40 x 1.25 ml	A325406
RealQ Plus 2x Master Mix Green Low ROX™	400 Reactions	4 x 1.25 ml	A324402
	4000 Reactions	40 x 1.25 ml	A324406
RealQ Plus 2x Master Mix Green Without ROX™	400 Reactions	4 x 1.25 ml	A323402
	4000 Reactions	40 x 1.25 ml	A323406