

## Nuclease-Free Water for molecular biology

MADE IN DENMARK

Nuclease-Free Water for molecular biology ID: 5404000	
Cat. No.	
A340037	1 x 100 ml
A340042	1 x 500 ml
A340044	1 x 1000 ml
A340099	Sample - 1 x 5 ml

### General Description

Nuclease-Free Water for molecular biology is suitable for use in PCR and other molecular biology techniques requiring nuclease-free conditions, such as processing of DNA or RNA (i.e. PCR, cDNA synthesis or qPCR).

- Water molecular biology grade
- Nuclease-free
- Tested free of DNase activity (endonuclease-, nicking and exonuclease activity)
- Tested free of RNase- and protease activity
- Absence of human DNA is tested by PCR
- Tested free of living organisms by a bioburden test

This product is highly purified and ultra-filtered. Not DEPC-treated.

Use this product as directed in the protocol of choice.

### Recommended Storage

Store at -20 °C to +25 °C.

Product expiry is stated on the label.

For Research Use Only. Not for use in diagnostics procedures.

Other product sizes, combinations and customized solutions are available. Please look at [www.ampliqon.com](http://www.ampliqon.com) or ask for our complete product list for PCR Enzymes. For customized solutions please contact us.

**Made in Denmark**

Issued 06/2025