

dATP – Single dNTP

Concentration: **100 mM**

MADE IN DENMARK

Cat. No.	dATP 100 mM ID: 5100910
A521102	1 x 250 µl

Key Features

- Purity \geq 99 % by HPLC.
- Ultrapure dNTPs for superior results in PCR and other applications

Ready-to-use molecular grade dNTP solution for use in DNA polymerization, DNA labeling and sequencing processes.

Recommended Storage and Stability

Long term storage at -20 °C. Product expiry at -20 °C is stated on the label.

Option: Store at +4 °C for up to 6 months.

Quality Control

The dNTPs are functionally tested with thermostable polymerases and tested for contaminating activities, with no traces of endonuclease activity, nicking activity and exonuclease activity.

Kit Components

dATP

100mM dATP in solution, pH 7.0 – 7.5

For Research Use Only. Not for use in diagnostics procedures.

Other product sizes, combinations and customized solutions are available. Please look at www.ampliqon.com or ask for our complete product list for PCR Enzymes. For customized solutions please contact us.

Made in Denmark

Issued 06/2025