

## dCTP – Single dNTP

Concentration: 100 mM

MADE IN DENMARK

Cat. No.	dCTP 100 mM ID: 5100920
A521202	1 x 250 µl

### Key Features

- Purity ≥ 99 % by HPLC.
- Ultrapure dNTPs for superior results in PCR and other applications

Ready-to-use molecular grade dNTP solution for use in DNA polymerization, DNA labeling and sequencing processes.

### Recommended Storage and Stability

Long term storage at -20 °C. Product expiry at -20 °C is stated on the label.

Option: Store at +4 °C for up to 6 months.

### Quality Control

The dNTPs are functionally tested with thermostable polymerases and tested for contaminating activities, with no traces of endonuclease activity, nicking activity and exonuclease activity.

### Kit Components

#### dCTP

100mM dCTP in solution, pH 7.0 – 7.5

For Research Use Only. Not for use in diagnostics procedures.

Other product sizes, combinations and customized solutions are available. Please look at [www.ampliqon.com](http://www.ampliqon.com) or ask for our complete product list for PCR Enzymes. For customized solutions please contact us.

**Made in Denmark**

Issued 06/2025