

## REAL-TIME PCR MASTER MIXES




### Real-time PCR Master Mixes with high efficiency and great prices


- All-in-one optimized real-time PCR Master Mixes.
- High efficiency for accurate quantification experiments.
- Convenient Hot Start capability for room temperature setup
- High specificity and stability.

RealQ Plus 2X Master Mixes from Ampliqon are easy to use all-in-one Master Mixes for both dye- and probe based detection. The mixes are available in three different levels of ROX™ reference dye and have been optimized by the Ampliqon R&D team to have high efficiency and specificity for more reliable results.


*Low Standard Deviation*



*High Efficiency*



*High Specificity*



	Size Reaction size: 25 µl	Content	Cat #
RealQ Plus 2x Master Mix Probe High ROX™	400 Reactions 4000 Reactions	4 x 1.25 ml 40 x 1.25 ml	A315402 A315406
RealQ Plus 2x Master Mix Probe Low ROX™	400 Reactions 4000 Reactions	4 x 1.25 ml 40 x 1.25 ml	A314402 A314406
RealQ Plus 2x Master Mix Probe Without ROX™	400 Reactions 4000 Reactions	4 x 1.25 ml 40 x 1.25 ml	A313402 A313406
RealQ Plus 2x Master Mix Green High ROX™	400 Reactions 4000 Reactions	4 x 1.25 ml 40 x 1.25 ml	A325402 A325406
RealQ Plus 2x Master Mix Green Low ROX™	400 Reactions 4000 Reactions	4 x 1.25 ml 40 x 1.25 ml	A324402 A324406
RealQ Plus 2x Master Mix Green Without ROX™	400 Reactions 4000 Reactions	4 x 1.25 ml 40 x 1.25 ml	A323402 A323406