

## IQON DNA LADDERS





### For direct gel loading and easy visualisation

- Cost saving solution
- Blue dye front running at 100-300 bp at 0.5 - 1.5 % agarose
- Suitable with TBE and TAE electrophoresis systems
- Clear and distinct bands provide easy estimation of PCR product lengths
- Provides easy size determination of double-stranded DNA

Three variants available:

1. **Iqon Mini DNA Ladder**, 100-500 bp, 100 lanes
2. **Iqon Low DNA Ladder**, 100-1000 bp, 100 lanes
3. **Iqon PCR Ladder**, 100-3000 bp, 100 lanes



	Size 5 µl/lane	Content	Cat #
Iqon Mini DNA Ladder, 100-500 bp	100 reactions	1 x 0.5 ml	A610441
Iqon Low DNA Ladder, 100-1000 bp	100 reactions	1 x 0.5 ml	A610541
Iqon PCR Ladder, 100-3000 bp	100 reactions	1 x 0.5 ml	A610641